

Absolute Diffusion Coefficients: Compilation of Reference Data for FCS Calibration

Peter Kapusta, PicoQuant GmbH, July 2010, Rev.1

Conventional FCS can provide absolute diffusion constant and concentration values only when the exact shape and size of the confocal volume is known. In practice, the shape is approximated by a 3D gaussian function. Size and elongation of the volume is obtained by calibration measurement of a reference compound with a known diffusion constant.

In order to get correct confocal volume parameters, the experimental conditions during calibration (wavelengths, excitation power, cover slip thickness, solvent and immersion medium) must be the same as during the measurement.

Precise diffusion coefficients of several reference compounds suitable for calibration purposes are summarized in the following table:

Fluorophore	λ_{Em} maximum in nm	Diffusion coefficient in water at 25°C (298.15 K) in $10^{-6} \text{ cm}^2\text{s}^{-1}$	Methods and references
Atto655-maleimid	686	4.07 ± 0.10 4.06 ± 0.09 4.09 ± 0.07	2fFCS [1], PFG-NMR [1] 2fFCS [4, 8] pmFCS [4]
Atto655-carboxylic acid	685	4.26 ± 0.08	2fFCS [1,3], PFG-NMR [1,3]
Atto655-NHS ester	685	4.25 ± 0.06	2fFCS [8]
Cy5	670	3.6 ± 0.1	2fFCS [8]
Alexa 647	665	3.3 ± 0.1	2fFCS [8]
Alexa 633	647	3.4 ± 0.1	2fFCS [8]
Rhodamine 6G	550	4.14 ± 0.05 4.3 ± 0.4 4.14 ± 0.01	2fFCS [1, 8] PFG-NMR [6] PB/CF [7]
Rhodamine B	560	4.5 ± 0.4 4.27 ± 0.04	PFG-NMR [6] PB/CF [7]
Rhodamine 123	530	4.6 ± 0.4	PFG-NMR [6]
Rhodamine 110	535	4.7 ± 0.4	PFG-NMR [6]
Fluorescein	520	4.25 ± 0.01	PB/CF [7]
Oregon Green 488	550	4.11 ± 0.06 4.10 ± 0.08	2fFCS [1] 2fFCS [8]
Atto488-carboxylic acid	523	4.0 ± 0.1	2fFCS [5]
TetraSpeck Beads, 0.1µm diameter	430 515 580 680	0.044 ± 0.07	2fFCS [2], DLS [2]

Abbreviations of measurement methods:

2fFCS, Dual Focus Fluorescence Correlation Spectroscopy; PFG-NMR, Pulsed Field Gradient Nuclear Magnetic Resonance; pmFCS, Polarization-Modulation Fluorescence Correlation Spectroscopy; PB/CF, Plug Broadening/Capillary Flow; DLS, Dynamic Light Scattering.

References:

- [1] Müller C. B., Loman A., Pacheco V., Koberling F., Willbold D., Richtering W., Enderlein J.: *Precise measurement of diffusion by multi-color dual-focus fluorescence correlation spectroscopy*. EPL, **83** (2008) 46001, <http://dx.doi.org/10.1209/0295-5075/83/46001>
- [2] Müller C. B., Weiß K., Richtering W., Loman A., Enderlein J.: *Calibrating differential interference contrast microscopy with dual-focus fluorescence correlation spectroscopy*. Optics Express, **16** (2008) 4322, <http://dx.doi.org/10.1364/OE.16.004322>
- [3] Dertinger T., Pacheco V., von der Hocht I., Hartmann R., Gregor I., Enderlein J.: *Two-focus fluorescence correlation spectroscopy: A new tool for accurate and absolute diffusion measurements*. ChemPhysChem, **8** (2007) 433, <http://dx.doi.org/10.1002/cphc.200600638>
- [4] Korlann Y., Dertinger T., Michalet X., Weiss S., Enderlein J.: *Measuring diffusion with polarization-modulation dual-focus fluorescence correlation spectroscopy*. Optics Express, **16** (2008) 4609, <http://dx.doi.org/10.1364/OE.16.014609>
- [5] Dertinger T., Ewers B.: *Unpublished results*. PicoQuant GmbH (2008)
- [6] Recalculated values from: Gendron P.-O. Avaltroni F., Wilkinson K. J.: *Diffusion coefficients of several rhodamine derivatives as determined by pulsed field gradient - nuclear magnetic resonance and fluorescence correlation spectroscopy*. Journal of Fluorescence, **18** (2008) 1093, <http://dx.doi.org/10.1007/s10895-008-0357-7>
- [7] Culbertson Ch. T., Jacobson S. C., Ramsey J. M.: *Diffusion coefficient measurements in microfluidic devices*. Talanta, **56** (2002) 365, [http://dx.doi.org/10.1016/S0039-9140\(01\)00602-6](http://dx.doi.org/10.1016/S0039-9140(01)00602-6)
- [8] Loman A., Müller C. B., Koberling F., Richtering W., Enderlein J.: *Absolute and precise measurements of the diffusion of small fluorescent dye molecules across the visible spectrum*. Poster, 14th International Workshop on Single Molecule Spectroscopy and Ultrasensitive Analysis in Life Sciences, September 17-19, 2008, Berlin, Germany

Note that the diffusion coefficient is not a constant. It is temperature dependent according to the Stokes-Einstein relationship,

$$D(T) = \frac{kT}{6\pi\eta(T)r}$$



where D is the diffusion coefficient, k is the Boltzmann constant ($1.3807 \cdot 10^{-23} \text{ J} \cdot \text{K}^{-1}$), r is the hydrodynamic radius. Moreover, the viscosity, η , of water or solvent in general is also a function of temperature. If the calibration experiment is not performed at $T=298.15 \text{ K}$, $t=25^\circ\text{C}$, the tabulated

values should be recalculated using this equation:

$$D(T) = D(25^\circ\text{C}) \cdot \frac{T}{298.15 \text{ K}} \cdot \frac{8.9 \cdot 10^{-4} \text{ Pa} \cdot \text{s}}{\eta(T)}$$

$$= D(25^\circ\text{C}) \cdot \frac{t + 273.15}{\eta(t)} \cdot 2.985 \cdot 10^{-6} \text{ Pa} \cdot \text{s} \cdot \text{K}^{-1}$$

The temperature dependence is strong and unfortunately often forgotten when comparing results. The effect is illustrated below for three popular fluorophores and the viscosity of water:



Helpful comments and suggestions of Prof. Wajih Al-Soufi, University of Santiago de Compostela, are gratefully acknowledged.

Copyright of this document belongs to PicoQuant GmbH. No parts of it may be reproduced, translated or transferred to third parties without written permission of PicoQuant GmbH. All Information given here is reliable to our best knowledge. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearances are subject to change without notice



PicoQuant GmbH
Rudower Chaussee 29 (IGZ)
12489 Berlin
Germany

Phone +49-(0)30-6392-6929
Fax +49-(0)30-6392-6561
Email info@picoquant.com
WWW www.picoquant.com