


FluoTime 200




PICOQUANT

High Performance Fluorescence Lifetime Spectrometer



- Modular and flexible design
- Picosecond time resolution
- Time-Correlated Single Photon Counting (TCSPC) via USB
- Laser injection optics
- Optional prism polarizers
- Advanced data analysis software



Applications

- Time-resolved fluorescence spectroscopy
- Fluorescence anisotropy decay analysis
- Ultra sensitive analysis
- Photochemistry
- Solar cell research
- Singlet oxygen
- Material research

Configuration and Standard Components

The FluoTime 200 spectrometer is a high performance fluorescence lifetime system featuring modular construction, single photon timing sensitivity and research flexibility. It contains the complete optics and electronics for recording fluorescence decays by means of Time-Correlated Single Photon Counting (TCSPC). The system can be used with femtosecond or picosecond lasers, like the picosecond diode lasers from PicoQuant. With the FluoTime 200, decay times down to a few picoseconds can be resolved. The system allows operation at laser repetition rates as high as 85 MHz and count rates up to several million counts/sec.



Laser coupling module

The laser coupling module consists of an iris, a beam steering mirror, an optional replaceable polarizer module and a focusing lens. The iris allows precise adjustment of the excitation intensity and prevents the penetration of room light into the sample chamber. The beam is directed towards the sample by an adjustable steering mirror. An optional prism polarizer selects the desired polarization plane of the non-polarized radiation. Finally, the beam is focused onto the sample by an adjustable lens.



Sample chamber

The sample chamber contains a versatile sample holder (standard: cuvette, optional: front face) and a mechanism to protect the detector from excessive light exposition. Temperature control of the cuvette holder is possible by attaching an external thermostat (tubing for the circulating fluid is pre-installed) or by replacing the standard cuvette holder with a peltier-cooled holder. Later upgrades to T- or X-geometry or motorized sample holders are possible.



Detector


The PMA detector unit, based on the Hamamatsu H5783 series photosensor module, is recommended for the majority of applications. The unit has a built-in high voltage power supply, signal pre-amplifier and a gold plated iron housing for optimal timing performance and maximum shielding. With this detector unit, an Instrument Response Function (IRF) as short as 200 ps can be achieved. Cooling is available for both detector types to reduce the number of dark counts.

Monochromator

For spectral filtering of the emitted light, a small monochromator is chosen. The standard single monochromator, with a 100 mm focal length, features low time dispersion, an $f/3.2$ aperture and externally adjustable slits with bandwidths of 4, 8 and 16 nm.

Data analysis


Numerical analysis of the raw data is an integral part of the TCSPC method. FluoFit is a Windows™ based global decay analysis software with an easy to use graphical user interface and presentation-ready numerical and graphical output. It implements an iterative reconvolution fitting routine with nonlinear error minimization. Various exponential decay models (up to fourth order) or rate constant distribution models can be fitted to the observed decay. FluoFit supports fitting of calculated anisotropy data as well. Running FluoFit simultaneously with the measurement software allows for parallel data analysis.



TCSPC data acquisition

The data acquisition module contains the complete timing electronics for Time-Correlated Single Photon Counting (TCSPC). The system works in forward start-stop mode, still operating at the full repetition rate of modern pulsed excitation sources. Several million photon events per second can be processed. Consequently, TCSPC acquisition times shorter than 1 second can be utilized. All functions of the TCSPC module are controlled by easy to use Windows™ based software. These functions include loading and saving experimental data as well as setting and storing measurement parameters. Optionally available is a data acquisition board for lower time resolution (phosphorescence decays).

Optical Layout



Options

Prism polarizers

For fluorescence anisotropy measurements a 15 mm aperture Glan-Laser polarizer can be provided. It features an excellent transmission range of 230-2700 nm and an extinction ratio of 10^{-5} .

Microchannel plate PMT

Detector based on the Hamamatsu R3809 series microchannel plate photomultiplier tube (MCP-PMT) to achieve the highest possible time resolution with the TCSPC method. With fast laser sources the IRF is less than 50 ps.

Multichannel scaler

For slow luminescence decays, the NanoHarp 250 multiscaler board can be added.

NIR-PMT

Detector based on the Hamamatsu H10330 series photomultiplier tube for fluorescence detection in the infrared up to 1400 nm (1700 nm as an option).

Versatile sample holder

The standard sample holder accommodates 1 cm (or smaller) path length cuvettes. Temperature control of the sample is made possible by circulating liquid and the tubing is pre-installed. An additional external thermostat (optional) can be attached. The cuvette holder can be optionally replaced with a peltier cooled cuvette holder or replaced with a simple front face sample holder.



Cryostat options



The Oxford Instruments *Optistat*TM series cryostat allows for low temperature photoluminescence measurements. Precise control of the sample temperature is possible with various cryostat types, ranging from 2.3 to 500 K.

Further options

Electrical stirrer, multi-pass sample holder, DLLs for custom programming, ...


Measurement Examples

Spectral dependence of the fluorescence decay kinetics (i.e. time-resolved fluorescence spectra) of an ethanol solution of Phenoxazon 9 (Nile Red) containing a trace amount of Coumarin 152. 32 decay curves were recorded (from 450 to 750 nm in 10 nm steps) using a FluoTime 200 with motorized monochromator. The decay curves show that, at longer times, the spectra are dominated by the Nile Red emission due to its longer lifetime. This measurement mode of the FluoTime 200 is particularly useful for studying solvation dynamics, excimer or exciplex formation and related phenomena.



Time-resolved fluorescence measurement of DASPI (10^{-5} molar in Ethanol) using a FluoTime 200 equipped with a MCP-PMT and the PicoHarp 300. The sample was excited with a LDH-P-C-405 picosecond diode laser with a 70 ps FWHM pulse width at 20 MHz repetition rate. Fluorescence response at 590 nm was detected through a single monochromator. The plot shows the measured Instrument Response Function (red), the sample decay (blue) and the fitted decay (black). The fitted fluorescence lifetime is 65 ± 5 ps.

The example on the right shows a time-resolved anisotropy measurement of Coumarin 153 (C153) and Coumarin 6 (C6) in dimethyl sulfoxide (DMSO) at a temperature of 290 K. The sample was excited with a LDH-P-C-405 diode laser and the emission at 500 nm was monitored with an MCP-PMT. The FluoTime 200 was equipped with Glan-Taylor prism polarizers. The figures show the recorded polarized decays, the IRF of the measurement and the calculated anisotropy decay together with a structure of the compounds. Fitting the anisotropy decay to a single exponential model without reconvolution yields a 110 ps and 330 ps rotational correlation time for C153 and C6, respectively. The results reflect the prolonged molecular shape of C6, in comparison to C153.



Specifications

System

Optical configuration L-Geometry (standard), T- or X-Geometry available
 Mode of operation Time-Correlated Single Photon Counting (TCSPC),
 Multichannel Scaling (MCS)¹⁾

Monochromators

Type single
 Focal length and aperture 100 mm F/3.2
 Stray light rejection $1 \cdot 10^{-5}$
 Gratings 1200 l/mm (standard)
 Step size (min) 0.5 nm
 Resolution 4 nm, 8 nm and 16 nm
 Wavelength accuracy 0.3 nm
 Wavelength repeatability 0.5 nm

Excitation Sources

Light source Pulsed LEDs (PLS Series) Laser Diode Heads (LDH Series) Ti:Sa Lasers
 Wavelengths 245-600 nm 375-510, 530, 635-1550 nm 240-950 nm²⁾
 Pulse width 400 ps - 1 ns 60-500 ps 20-200 fs
 Repetition rate up to 40 MHz up to 40 MHz (optional 80 MHz) up to 84 MHz

Detectors

Type³⁾ PMT (PMA Series) MCP-PMT NIR-PMT
 Spectral range 185-700 nm, 185-820 nm, 300-900 nm 160-650 nm, 160-850 nm 950-1400 nm
 Dark counts* (uncooled) < 50 cps < 900 cps < 6000 cps < 50 cps < 500 cps ---
 Dark counts* (cooled) --- < 50 cps < 900 cps < 10,000 cps

Data Acquisition Unit

Type PicoHarp 300 NanoHarp 250
 Time resolution (bin width) 4 ps 4 ns, 32 ns
 Count depth 16 bit 18 bit
 Dead time < 95 ns 8 ns
 Lifetime resolution < 10 ps < 10 ns
 Time channels per curve up to 65536 1024 to 262144
 Differential non-linearity < 1 % rms, < 5 % pp < 0.1 % rms, < 1 % pp.
 Discriminator level, zero cross adjust software adjustable
 Collection time 1 ms - 10 h

Data Analysis Software

Type FluoFit Pro
 Analysis possibilities exponential decay, lifetime distribution, anisotropy, global analysis

Operation

Operating system Windows™ 2000/XP/Vista/7

Electrical and Dimensional

Power requirements 220/240 or 110/120 VAC, 50/60 Hz
 Dimensions (standard) 800 × 400 × 300 mm (w × d × h)

¹⁾ with a different photon counting board, ²⁾ incl. harmonics, ³⁾ other detectors and cooling available upon request, *typ. value at 20 °C ambient temperature

Other PicoQuant Systems

FluoTime 300 "EasyTau"
 Automated fluorescence
 lifetime spectrometer




MicroTime 100
 Upright time-
 resolved
 fluorescence
 microscope



MicroTime 200
 Inverse time-resolved
 fluorescence
 microscope



LSM Upgrade Kit
 Compact lifetime & FCS upgrade
 kit for Laser Scanning Microscopes



All Information given here is reliable to our best knowledge. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearances are subject to change without notice. Trademarks or corporate names are used for explanation and identification, to the owner's benefit and without intent to infringe.

© PicoQuant GmbH, December 2012



PicoQuant GmbH
 Rudower Chaussee 29 (IGZ)
 12489 Berlin
 Germany

Phone +49-(0)30-6392-6929
 Telefax +49-(0)30-6392-6561
 Email info@picoquant.com
 WWW http://www.picoquant.com