

# Luminosa

## Single Photon Counting Confocal Microscope


- Software-controlled confocal system based on an inverted microscope
- Versatile excitation system with laser wavelengths from 375 to 1064 nm
- VarPSF: fine-tuned observation volume for FCS and single molecule FRET experiments
- Motorized positioning stage for “tiling and stitching” in transmission, epi, DIC and confocal modes
- Scanning options: FLIMBee galvo scanner and piezo objective scanning
- Stage top incubator for extended live cell experiments
- Up to 6 truly parallel detection channels with SPAD and/or Hybrid PMTs
- PDA-23 add-on: an integrated SPAD array detector
- < 650 ps dead time per channel and 5 ps time bins
- One click auto-alignment for consistent optimal performance
- Fast results with minimal user interaction thanks to GPU accelerated algorithms and context-based workflows for FCS, FLIM, and single molecule detection
- Imaging scanning microscopy (ISM-FLIM): FLIM with more spatial resolution and higher contrast
- An open Python interface enabling adaptive workflows, custom acquisition and reporting

### Core Methodologies

- Single molecule FRET (burst and time trace analysis)
- Fluorescence lifetime imaging (FLIM)
- Fluorescence correlation spectroscopy (FCS)
- FLIM-FRET
- Anisotropy imaging
- Fluorescence recovery after photobleaching (FRAP)
- Time-resolved photoluminescence
- Antibunching

### Application Areas

- Dynamic structural biology
- Studying cellular mechanisms driven by phase separation
- Environmental sensing and marker multiplexing
- Mapping dynamics and structure of cellular membranes
- Characterizing functional nanovesicles
- Metabolic Imaging
- Studying chemical reactions at the single molecule level
- DNA nanotechnology
- Characterizing advanced materials



## Specifications

Feature Highlights	
LumiFinder	<ul style="list-style-type: none"> <li>Automated detection and measurement of immobilised emitters</li> <li>Automatically use different excitation power for emitter identification and time trace measurements to avoid unnecessary photobleaching</li> <li>Advanced stopping criteria for the time trace measurement offer more flexibility in single molecule FRET measurements of immobilised emitters</li> </ul>
FRETcompass	Automated determination of correction factors used in single molecule FRET experiments based on the benchmark study published in 2018 (Hellenkamp et al. Precision and accuracy of single-molecule FRET measurements – a multi-laboratory benchmark study. Nat Methods 15, 669-676 (2018))
InstaFCS	Online FCS fitting with automatic suggestions for model and parameters
InstaFLIM	<ul style="list-style-type: none"> <li>Provides first suggestions for FLIM species separation without user interaction</li> <li>Suggested parameters can be fine-tuned</li> <li>Simultaneous TCSPC and phasor analysis options for ROI determination</li> </ul>
Sample-free Auto-alignment (SFA)	<ul style="list-style-type: none"> <li>One-click auto-alignment without any sample required</li> <li>Ensures optimal performance for every measurement, even the most challenging ones</li> </ul>
CalEx*	<ul style="list-style-type: none"> <li>Excitation laser power calibration allows setting and displaying excitation intensity in <math>\mu\text{W}</math></li> <li>No external power meter required</li> </ul>
VarPSF*	Switch between diffraction limited and larger observation volume to: <ul style="list-style-type: none"> <li>Increase observation time window to check for dynamic transitions between states</li> <li>Study diffusion properties of larger particles, 3-6 times larger volumes are achieved (depending on wavelength)</li> </ul>
LumiPy	Enabling adaptive workflows and custom acquisition and reporting via a python interface allowing for Luminosa functionality to be exported to Python classes

Excitation	
	1 - 8 lasers with choice of wavelengths from 375 - 1064 nm <ul style="list-style-type: none"> <li>Optional: 2 additional exit ports available for coupling additional lasers</li> <li>Optional: fs fiber laser from TOPTICA at 780 nm for metabolic imaging</li> </ul>
Commonly used wavelengths	405, 440, 485, 532, 560, 595 and 640 nm
Individual laser operation modes	Pulsed and Continuous Wave (continuous mode availability depends on the laser model)
Repetition rates	1 - 40 MHz, for specific wavelengths up to 80 MHz
Laser driver	8 channel Sepia PDL 828 L
Supported multi-line operation	Pulsed Interleaved Excitation (PIE) or simultaneous emission modes (for specific wavelengths alternating excitation move after (ALEX) in $\mu\text{s}$ and $\text{ms}$ time scales also available)
Motorized main dichroic wheel	6 positions available for 18 × 25 × 3 mm single, dual, triple, and quad-band dichroics

Microscope	
	Inverted Evident IX73 (partially motorized)
Objective types	<ul style="list-style-type: none"> <li>UPLXAPO60XW X-Line PlanApochromat, NA 1.2, water immersion, 400 - 900 nm</li> <li>UPLXAPO100x1.45 Oil X-Line PlanApochromat, NA 1.4, oil immersion, 400 - 850 nm</li> <li>Other oil-immersion, apochromatically corrected, air-spaced, NIR-enhanced, or long-working-distance objectives are available</li> </ul>
Observation modes	<ul style="list-style-type: none"> <li>Confocal, transmission</li> <li>Optional: epi-widefield illumination, differential Interferometric Contrast (DIC)</li> </ul>
Oculars	Yes
Camera	Optional CCD camera at the left side port
Sample holder	<ul style="list-style-type: none"> <li>Multi-functional sample holder compatible with several standard sizes of microscope coverslips and petri dishes</li> <li>96-well plate sample holder</li> </ul>
Stage-top incubator	Yes
Axial drift (typical values **)	<ul style="list-style-type: none"> <li>For 15 minutes measurement: <math>50 \pm 10</math> nm</li> <li>For 2 hour measurement: <math>100 \pm 20</math> nm</li> </ul>

Detection Options		
Point detection options	<ul style="list-style-type: none"> <li>Parallel detection with up to 6 point detectors in the spectral range of 300 - 1000 nm</li> <li>23 pixel SPAD array</li> <li>Additional fiber coupled or free-space exit ports for customised detectors</li> </ul>	
Point detection types	<ul style="list-style-type: none"> <li>Single Photon Avalanche Diodes (SPAD) for low dark-counts and maximum detection efficiency</li> <li>Hybrid photomultiplier tubes (PMA-Hybrid-40) for optimum timing performance and when working at high count rates</li> </ul>	
Type	SPAD (SPCM-AQRH)	PMA Hybrid - 40
Spectral range	400 - 1000 nm	300 - 720 nm
Dark counts (at 20 °C, typ. value)	< 250 cps	< 300 cps
Photon detection efficiency (PDE)	Typ. 50 % at 550 nm Typ. 70 % at 650 nm	Typ. 40 % at 500 nm
Timing response	200 - 400 ps (depending on wavelength)	< 120 ps
PDA-23 Array Detector with 23 independent Single Photon Avalanche Diode pixels		
Active area	100 $\mu$ m x 110 $\mu$ m	
Effective fill factor with microlenses	> 75 %, typical 80 %	
Photon detection efficiency (PDE)	<ul style="list-style-type: none"> <li>Typ. peak of PDE 55 % at 520 nm</li> <li>PDE typically above 35 % for wavelengths from 440 to 660 nm</li> <li>Usable spectral window 400 - 950 nm</li> </ul>	
Dark counts at working temperature of 10 °C via integrated cooling	<ul style="list-style-type: none"> <li>Typ. 100 cps for most pixels</li> <li>Median value of dark counts per second is &lt; 100</li> <li>Typ. 1 hot pixel with &gt; 1 kcps and below &lt; 10 kcps</li> <li>Hot pixel not in the center of the array</li> </ul>	
Jitter	< 120 ps FWHM for wavelengths between 510 and 785 nm	

Deadtime of each pixel	50 ns
Crosstalk probability between pixels	0.1 %
Afterpulsing probability	0.1 %
Maximum count rate per pixel	Typically > 7 Mcounts per second

Data Acquisition	
	<ul style="list-style-type: none"> <li>Based on Time-Correlated Single Photon Counting (TCSPC) in Time-Tagged Time Resolved (TTTR) measurement mode in T3 and T2 modes</li> <li>Simultaneous data acquisition with up to 6 channels + PDA-23</li> <li>New dynamic binning format for FLIM imaging</li> </ul>
TCSPC units	<ul style="list-style-type: none"> <li>MultiHarp 150 P (standard option)</li> <li>MultiHarp 160 (for PDA-23)</li> </ul>
Time resolution (bin width)	5 ps
Dead time	< 650 ps
Timing precision	< 28 ps RMS
Number of available detector inputs	8 for MultiHarp 150 P (standard option) 32 for MultiHarp 160 (for PDA-23)

Scanning Options	
FLIMbee XY galvo scanner	<ul style="list-style-type: none"> <li>Min. pixel dwell time: 0.5 <math>\mu</math>s,</li> <li>Min. pixel size: 17 nm (60x objective)</li> <li>Up to 5.2 frames per second for a 512 x 512 pixel image max.</li> <li>Scan field: 200 x 200 <math>\mu</math>m (60x objective)</li> <li>Possibility to bypass the FLIMbee galvo scanner for point measurements</li> </ul>
XY objective scanning	<ul style="list-style-type: none"> <li>Min. pixel dwell time: 0.2 ms min.</li> <li>Pixel size: 1 nm (independent of objective)</li> <li>Max. scan field: 80 x 80 <math>\mu</math>m (independent of objective)</li> </ul>
XY-scanning combined option	Combination of FLIMbee and objective scanning on the same system with software-controlled change between the two scanning modalities
Piezo-based Z-scanning	<ul style="list-style-type: none"> <li>Min. step size: 50 nm</li> <li>Overall range: 100 <math>\mu</math>m (optional 400 <math>\mu</math>m)</li> </ul>
Positioning stage with spiral scan mode, creating overview map by tiling and stitching	<ul style="list-style-type: none"> <li>Overall range: 121 x 81 mm</li> <li>Max speed: 300 mm/second</li> <li>Positioning repeatability: &lt; 0.15 <math>\mu</math>m</li> </ul>

Luminosa Software	
Operating System	Windows 11
GPU-based programming	Software-Open GL
Context based workflows for easy acquisition and analysis	For single molecule FRET burst analysis, single molecule imaging, FCS in vitro, FCS in cells, FLIM, FLIM-FRET, steady-state FRET, anisotropy imaging, and FRAP

General Info	
Operating environment	<ul style="list-style-type: none"> <li>Room temperature range 15 - 25 °C</li> <li>Room temperature stability <math>\pm 1.5</math> °C (recommended)</li> <li>Room humidity &lt; 60 %</li> </ul>
Optical table (minimum dimensions)	<ul style="list-style-type: none"> <li>For 1 - 4 detectors: 1500 × 900 mm</li> <li>For 5 - 6 detectors: 2000 × 900 mm</li> </ul>
Power consumption	<ul style="list-style-type: none"> <li>6 A at 230 V AC (typ. EU)</li> <li>20 A at 110 V AC (typ. USA)</li> </ul>
Operating voltage	115 or 230 V AC
Altitude	Guaranteed performance up to 2000 m above sea level

\* Feature available only for PicoQuant lasers

\*\* Under constant room temperature and after the sample and z drive have settled for at least an hour



PicoQuant GmbH  
 Rudower Chaussee 29 (IGZ)  
 12489 Berlin  
 Germany

Phone +49-(0)30-1208820-0  
 Email [info@picoquant.com](mailto:info@picoquant.com)  
 Web [www.picoquant.com](http://www.picoquant.com)