



LDH-P-FA-530



PICOQUANT
Unternehmen für optoelektronische
Forschung und Entwicklung

<http://www.picoquant.com>

Picosecond Pulsed Laser Diode Heads at 530 nm



- Center wavelength: 531 ± 3 nm
- Pulse width below 100 ps (FWHM)
- Up to 2 mW average output power
- Repetition rate from 10 kHz to 40 MHz
- FC/APC fibre connector



Applications

- Time-resolved fluorescence spectroscopy & microscopy (FLIM, FRET, FCS)
- Biochemical analytics
- Time-response characterization of opto-electronic devices
- Quantum cryptography
- Diffuse Optical Tomography (DOT)
- Fibre characterization

