

collimating lens (standard)

10 mm VIS
15 mm UV

filter (optional)




10 mm color glass filter
15 mm band pass filter

10/15

10/15

51

approx. 1400
cable



PLS series		Dimension sheet		General Dimensions without tolerance indication		all dimensions [mm]		Scale: M 1:1					
Drawing Standard DIN 406				<table border="1"> <tr> <th>Name</th> <th>Date</th> </tr> <tr> <td>Kulisch</td> <td>05.10.2018</td> </tr> </table>		Name	Date	Kulisch	05.10.2018	Last Update: 05.10.2018		Weight: 166.34 g	
Name	Date												
Kulisch	05.10.2018												
				The content of this drawing is the intellectual property of Picoquant and may not be copied or forwarded to third parties									
						Title:							
				PICOQUANT Rudower Chaussee 29 D-12489 Berlin Tel: +49-30-1208820-0 Fax: +49-30-1208820-90 info@picoquant.com		Drawing Number PQ.026.1010.0001.12		Sheet 4 / 4					
						Supplier No:		Feedstock:					
								Size: A3					