

# Mapping Molecules Quantitatively in Confocal Fluorescence Microscopy

Marcelle König, Caroline Berlage, Paja Reisch, Felix Koberling, Haisen Ta,  
Rainer Erdmann

PicoQuant GmbH, Rudower Chaussee 29, 12489, Berlin, Germany



PICOQUANT

24th International Workshop on  
“Single Molecule Spectroscopy and Super-resolution Microscopy in the Life Sciences”  
September 12-14, 2018 in Berlin, Germany

# Quantification of Molecules

**REPORT**

**Counting Low-Copy Number Proteins in a Single Cell**

Bo Huang<sup>1\*</sup>, Hongkai Wu<sup>1,†</sup>, Devaki Bhaya<sup>2</sup>, Arthur Grossman<sup>2</sup>, Sebastien Granier<sup>2</sup>, Brian K. Kobilka<sup>3</sup>, Richard N. Zare<sup>1,‡</sup>

<sup>\*</sup> Present address: Department of Chemistry and Chemical Biology, Harvard University, Cambridge, MA 02138, USA.

<sup>†</sup> Present address: Department of Chemistry, Tsinghua University, Beijing 100084, China.

<sup>‡</sup> See all authors and affiliations

Science 05 Jan 2007  
Vol. 315, Issue 5806, pp. 81-84  
DOI: 10.1126/science.1133992

Editorial

**NATURE METHODS | VOL.9 NO.7 | JULY 2012**

The quest for quantitative microscopy

With the aid of informatics, microscopy is in the evolution into a more quantitative and powerful

Stochastic protein expression in individual cells at the single molecule level

Long Cai, Nir Friedman & X. Sunney Xie

362 (16 March 2006)

Received: 12 September 2005  
Accepted: 23 January 2006  
Published online: 16 March 2006


**MB&C | TECHNICAL PERSPECTIVE**

**Every laboratory with a fluorescence microscope should consider counting molecules**


Valerie C. Coffman<sup>a</sup> and Jian-Qiu Wu<sup>a,b</sup>

<sup>a</sup>Department of Molecular Genetics and <sup>b</sup>Department of Molecular and Cellular Biochemistry, The Ohio State University, Columbus, OH 43210

→ Variety of research objectives for counting and mapping of molecules and their concentrations in cells



- Limited to small numbers
- Destructive



- Needs calibration measurements

### Fluctuation Analysis Based Methods


- FCS / FLCS
- PCH
- N&B
- ...

### Localization Microscopy

- PALM
- dSTORM
- PAINT
- ...

# Outline

## THE METHOD: Counting by Photon Statistics (CoPS)



## PROOF OF PRINCIPLE: Measurements with Origami



### ARTICLE

Received 24 Mar 2014 | Accepted 2 Jul 2015 | Published 13 Aug 2015



DOI: 10.1038/ncomms8977

OPEN

## Mapping molecules in scanning far-field fluorescence nanoscopy

Haisen Ta<sup>1</sup>, Jan Keller<sup>1</sup>, Markus Haltmeier<sup>2,3</sup>, Sinem K. Saka<sup>4</sup>, Jürgen Schmied<sup>5</sup>, Felipe Opazo<sup>4</sup>, Philip Tinnefeld<sup>5</sup>, Axel Munk<sup>2,6</sup> & Stefan W. Hell<sup>1</sup>

## TOWARDS BIOLOGICAL SAMPLES: First results




# Counting by Photon Statistics (CoPS)

The Principle behind  
**Counting by Photon Statistics (CoPS)**  
 is similar to antibunching:


**A single molecule can only emit one photon at a time.**

Method developed by Dirk-Peter Herten,  
 Heidelberg University

**Confocal microscope with pulsed excitation and four detectors**




**Detection of coincident photons**  
 (photons that arrive after the same laser pulse)



Adapted from Grußmayer et al.,  
 Phys. Chem. Chem. Phys., 2017, Suppl.

**Measurement of the distribution of multiple photon detection events**

Relative probabilities depend on number of emitters  $N$ , individual brightness  $p$  and number of detectors  $m$ .




Adapted from Grußmayer et al.,  
 Phys. Chem. Chem. Phys., 2017, Suppl.


# Mapping Molecules based on Counting by Photon Statistics (CoPS)

## Multi-photon detection events (immobilized DNA Origami with 9 ATTO647N)

1-photon events




2-photon events




scale bar: 1  $\mu\text{m}$


3-photon events




4-photon events



Estimated molecular brightness  $p_{\text{est}}$




Estimated emitter density  $n_{\text{est}}$



Method published by Haisen Ta *et al.*, Nature Communications, 2015

# Map of Molecule Distributions

## Analysis of multi-photon detection events (immobilized Origami with 9 ATTO647N)




Grey scalebar: Intensity [photons /pixel];  
Color scalebar: Number of emitters per spot (summed up density per pixel)

Parameters:

10  $\mu$ W excitation, 300  $\mu$ s px dwell time, 20nm px size, 10 MHz, 500 x 500 px


## Histogram of the number of emitters in one origami



# Proof of Principle: Red DNA Origami

## Red DNA-Origami with varying number of emitters (GattaQuant)

- 1 ATTO647N
- 4 ATTO647N
- 9 ATTO647N
- 17 ATTO647N
- 23 ATTO647N
- 30 ATTO647N




[http://www.gattaquant.com/files/gatta-brightness\\_product\\_sheet\\_1.pdf](http://www.gattaquant.com/files/gatta-brightness_product_sheet_1.pdf)

## Expected numbers of emitters per origami:

Calculation assuming binomial distribution with


- $n$  binding sites
- binding probability  $p$



# Proof of Principle: Red DNA Origami

## Red DNA-Origami with varying number of emitters

- 1 ATTO647N
- 4 ATTO647N
- 9 ATTO647N
- 17 ATTO647N
- 23 ATTO647N
- 30 ATTO647N




[http://www.gattaquant.com/files/gatta-brightness\\_product\\_sheet\\_1.pdf](http://www.gattaquant.com/files/gatta-brightness_product_sheet_1.pdf)

## Expected numbers of emitters per origami:

Calculation assuming binomial distribution


## Measured brightness for increasing numbers of emitters per origami:

- Number of detected photons per identified origami in image
- Normalized for one emitter




# Counting by Photon Statistics: Results with Red DNA Origami


## 1 Emitter




## 4 Emitters




## 9 Emitters




## 17 Emitters




## 23 Emitters



## 30 Emitters



# Red DNA Origami: Overestimation for Higher Numbers



CoPS **overestimates** the emitter number for higher numbers per cluster.


Possible issues:

- **Saturation** of detection electronics
- Detector **afterpulsing**
- **Interaction** of fluorophores in DNA origami

# Proof of Principle: Blue/Green DNA Origami

## Blue/green DNA-Origami with varying number of emitters (GattaQuant)

- 1 ATTO488
- 4 ATTO488
- 12 ATTO488
- 24 ATTO488




### Expected numbers of emitters per origami:

Calculation assuming binomial distribution with

- $n$  binding sites
- binding probability  $p$


### Measured brightness for increasing numbers of emitters per origami:

- Number of detected photons per identified origami in image
- Normalized for one emitter




# Counting by Photon Statistics: Results with Blue/Green DNA Origami


## 1 Emitter




## 4 Emitters



## 12 Emitters




## 24 Emitters





# Biological Samples


## Plasma membranes with YFP S. Munck, KU Leuven



## Additional challenges in biological samples


- Even higher background
- Very dense, overlapping clusters
- Not two-dimensional, differences in z-position of clusters
- Fluorescent proteins not as bright, but slightly more stable than Atto488

## Nuclear Pore Complex with eGFP A. Rybina, A. Politi, J. Ellenberg, EMBL, Heidelberg




# Biological Samples: Nuclear Pore Complex (16 Emitters expected, EGFP)

Bottom of single interphase cell: Single pores with 16 emitters each




Calculated emitter density per pixel



scale bar: 1  $\mu$ m


Histogram of emitter numbers per pore




Homozygous cell line NUP214-mEGFP  
(Sample kindly provided by Arina Rybina, Antonio Politi, Jan Ellenberg, EMBL)

# Summary of What Works So Far...

## Data Acquisition with MicroTime200, four SPADs and HydraHarp



## Analysis with Matlab-Software from Haisen Ta



## Sample Requirements

- Fixed sample with low background
- Bright and stable fluorophores (preferably 640 nm excitation)
- Quantification of single clusters
- Less than 10 Emitters
- Narrow distribution of emitter numbers
- 2D

Interested?


Please contact us!



PICOQUANT

# Acknowledgement


## Haisen Ta




**Dirk-Peter Herten** and lab members  
**Johan Hummert** and **Wioleta Chmielewicz**, Heidelberg University

## Sample preparation

• **Sebastian Munck**  
KU Leuven



• **Arina Rybina,**  
**Antonio Politi,**  
**Jan Ellenberg,**  
EMBL, Heidelberg



**Caroline Berlage**, M.Sc. student  
(Supervisor **Oliver Benson**, HUB)

Poster P2: Molecular Counting  
by Photon Statistics

- Experimental parameters
- Origami with Atto488
- Photobleaching
- Limitations and Outlook