

collimating lens (standard)

10 mm VIS
15 mm UV

filter (optional)




10 mm color glass filter
15 mm band pass filter

10/15

10/15

51

approx. 1400
cable



PLS series		Dimension sheet		General Dimensions without tolerance indication		all dimensions [mm]		Scale: M 1:1	
Drawing Standard DIN 406				Date	Last Update: 05.10.2018			Weight: 166.34 g	
				05.10.2018					
The content of this drawing is the intellectual property of Picoquant and may not be copied or forwarded to third parties									
PICO QUANT Rudower Chaussee 29 D-12489 Berlin Tel: +49-30-1208820-0 info@picoquant.com						Title:			
						Drawing Number PQ.026.1010.0001.12		Sheet 4 / 4	
						Supplier No:		Feedstock:	
								Size: A3	