



**WARNING - LASER RADIATION**  
 AVOID EXPOSURE TO BEAM  
 CLASS 3B LASER PRODUCT  
 Maximum Power  $P_o < 100 \text{ mW}$   
 Wavelengths  $\lambda = 450 \text{ nm}$   
 SEE MANUAL  
 Complies with IEC 60825-1:2014/21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007.

PicoQuant GmbH  
 Rudower Chaussee 29  
 D-12489 Berlin  
 Germany  
 Model No. Unico-450  
 Part No. 917450  
 Serial No.  
 Mfg. Date:  
 www.picoquant.com

Unico Drawing Standard DIN 406 	Dimension sheet	General Dimensions without tolerance indication		Scale: M 1:1
		Date 17.12.2024	Weight: approx. 953.64 g	
The content of this drawing is the intellectual property of Picoquant and may not be copied or forwarded to third parties		Title: PQ.088.1000.1101.12_Unico_Webpage		
PICO QUANT Rudower Chaussee 29 D-12489 Berlin Tel: +49-30-1208820-0 info@picoquant.com		Drawing Number PQ.088.1000.1001.12	Supplier No: N.N. Feedstock: N.N.	Sheet 1 / 1 Size: A3